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
Uras Techno Vibro Pot® permits high-speed grinding of small-quantity powder of any kind, from grain to ultra fine powder. Especially, as a grinding mill for fine ceramics etc., various drum materials are provided: stainless steel, porcelain and natural rubber-lining.

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
●Principle of Operation

Vibro Pot® is a vibration utilizing test grinding mill for efficient grinding and mixing. It makes use of forced circular vibration to impart motion to a drum filled with material to be ground and grinding media (balls or rods).

Therefore, the material is ground by impact, shearing and friction action with the media and the drum wall. The grinding effect, having vibratory acceleration of $50\sim 100\text{m/s}^2$ (5~10G), shortens the necessary time to less than 1/10 of the time compared to that of rotating drum type mill.

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- Standard voltage is 200V class, but 400V class are also available.
- Able to operate at vibrating acceleration 100m/s^2 (10G), if media size is less than φ12mm.
- All kind of media balls mentioned following.
If you use rods as media please contact us.

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
● Features

1. With strong vibratory acceleration up to 100m/s² (10G), it can produce a remarkable grinding effect and is suitable for grinding material down to ultra fine powder of 1 μm or less.
2. Compact size and easy maintenance are realized by the use of Uras Vibrator of an established reputation.
3. Dry and Wet Grinding
Both dry and wet grinding are possible. A wide selection of drum materials are: stainless steel, porcelain and natural rubber-lining.
4. The standard construction employing a soundproofing cover causes negligible environmental disturbance.
5. Mixing Function
Besides grinding, instant homogeneous mixing is possible.

● Application

- Metal compounds, such as ferrite, barium titanate, molybdenum disulfide and alumina.
- Fire proof materials, such as silica stone, chamotte and bauxite.
- Food, medicine and organic materials.
- Fine ceramics and other new materials.

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形 式 Type	寸 法 Dimension mm										質 量 Mass kg
	A	H	I	J	K	L ₁	L ₂	M	N	P	
YAMP-2SND	340	600	460	100	40	740	450	700	560	440	120
-4JND											130
-6SND	350	640	500	100	40	840	600	800	710	590	175
-6JND											185
-6GND											185

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